



Now with Isolator Wrap™

Environmental Monitoring Media Confidence. It's Critical.



- Superior Quality
- Unsurpassed Sterility Assurance
- Excellent Service

remel

Part of Thermo Fisher Scientific

Providing the Critical Confidence you need in all your

Now with Isolator Wrap™

- Foil barrier wrap impermeable to isolator decontamination process
- An extra bag for transport to the laboratory
- Irradiated GripFit™ plates with labeling on the bottom of plate
- Extended shelf life - 6 months from date of manufacture

Superior Quality

- The highest quality raw materials and components are used in production
- Each lot is tested to ensure consistent, reproducible growth performance
- Large, single lots available to minimize time spent on QC testing in your lab

Exceptional Sterility Assurance

- Manufactured in a Class 10,000 (ISO Class 7) cleanroom with Class 100 (ISO Class 5) work zones to ensure high levels of sterility
- Sterile plated media receives gamma-irradiation treatment validated in accordance with ISO 11137-1,2,3 standards to achieve the highest available sterility assurance level (SAL= 10^{-6})

Excellent Service

- Reliable, next-day delivery provides the media you need when you need it
- Dedicated Customer Service Representatives are ready to help with your order
- Expert Technical Support available to answer your specific questions
- Lot-specific Certificates of Quality are available online 24/7 for your convenience



ur controlled environments.



STERILE CONTACT PLATES

Standard Sterile Contact Plates and Isolator Wrap™ Sterile Contact Plates are manufactured in a cleanroom environment and gamma-irradiated to ensure sterility [SAL=10⁻⁶]

STANDARD CONTACT PLATES, STERILE* are double bagged with a third bag enclosed for customer use (Double + Bag).

DESCRIPTION	TYPE/FILL	UNIT	REF
Contact Plate, Sterile, D/E Neutralizing Agar (Irradiated)	Double + Bag	10/pk	R111803
	Double + Bag	100/pk	R111804
Contact Plate, Sterile, Sabouraud Dextrose Agar with Lecithin, Polysorbate 80 (pH 5.6) (Irradiated).	Double + Bag	10/pk	R111805
	Double + Bag	100/pk	R111806
Contact Plate, Sterile, Tryptic Soy Agar with Lecithin, Polysorbate 80 (Irradiated)	Double + Bag	10/pk	R111800
	Double + Bag	100/pk	R111802
Contact Plate, Sterile, Custom†	Double + Bag	Made to Order	

ISOLATOR WRAP™ STERILE CONTACT PLATES* are contained in a protective outer barrier wrap. An extra bag is provided for transport of plates to the laboratory for testing. For use in isolator environments.

DESCRIPTION	TYPE/FILL	UNIT	REF
Isolator Wrap™ Sterile D/E Neutralizing Agar	Foil barrier wrap + bag	10/pk	R111823
	Foil barrier wrap + bag	100/pk	R111824
Isolator Wrap Sterile Sabouraud Dextrose Agar with Lecithin, Polysorbate 80 (pH5.6)	Foil barrier wrap + bag	10/pk	R111825
	Foil barrier wrap + bag	100/pk	R111826
Isolator Wrap Sterile Tryptic Soy Agar with Lecithin, Polysorbate 80	Foil barrier wrap + bag	10/pk	R111820
	Foil barrier wrap + bag	100/pk	R111822

*Not intended for IVD use. † GripFit™ not available on these products.



CONTACT I* (Single Bag) and CONTACT II* (Double Bag) are manufactured in a cleanroom environment.

DESCRIPTION	TYPE/FILL	UNIT	REF
Contact I Blood Agar (TSA w/5% Sheep Blood)	Single Bag	20/pk	R111007
Contact I D/E Neutralizing Agar	Single Bag	20/pk	R111103
	Single Bag	100/pk	R111104
Contact I Lethen Agar [†]	Single Bag	20/pk	R111019
Contact I Malt Extract Agar	Single Bag	20/pk	R111021
Contact I Mycological Agar w/Lecithin, Polysorbate 80 [†] (Mycophil Agar)	Single Bag	20/pk	R111030
Contact I Potato Dextrose Agar [†]	Single Bag	20/pk	R111045
Contact I Rose Bengal Agar w/chloramphenicol	Single Bag	20/pk	R111042
	Single Bag	100/pk	R111105
Contact I Sabouraud Dextrose Agar w/Lecithin, Polysorbate 80 (pH 5.6)	Single Bag	20/pk	R111105
	Single Bag	100/pk	R111106
Contact I Sabouraud Agar, Emmons	Single Bag	20/pk	R111051
Contact I Tryptic Soy Agar	Single Bag	20/pk	R111083
Contact I Tryptic Soy Agar w/Lecithin, Polysorbate 80	Single Bag	20/pk	R111101
	Single Bag	100/pk	R111102
Contact II D/E Neutralizing Agar	Double Bag	20/pk	R111203
	Double Bag	100/pk	R111204
Contact II Sabouraud Dextrose Agar w/Lecithin, Polysorbate 80 (pH 5.6)	Double bag	20/pk	R111205
	Double bag	100/pk	R111206
Contact II Tryptic Soy Agar w/Lecithin, Polysorbate 80	Double bag	20/pk	R111201
	Double Bag	100/pk	R111202

SETTLING PLATES, STERILE* are double bagged with a third bag enclosed for customer use (Double + Bag). They are manufactured in a cleanroom and gamma-irradiated for sterility (SAL=10⁻⁶).

DESCRIPTION	TYPE/FILL	UNIT	REF
Settling Plate, Sterile, Tryptic Soy Agar (Irradiated)	Monoplate - Double + Bag	10/pk	R111870
	Monoplate - Double + Bag	100/pk	R111872
	15x150mm - Double + Bag	20/pk	R111875
Settling Plate, Sterile, Custom	Double + Bag	Made to Order	

SETTLING PLATES* (Single and Double Bag) are manufactured in a cleanroom environment.

DESCRIPTION	TYPE/FILL	UNIT	REF
Settling Plate, Tryptic Soy Agar	Monoplate - Single Bag	10/pk	R111570
	15x150mm - Single Bag	10/pk	R111575
	Monoplate - Double Bag	10/pk	R111670
	Monoplate - Double Bag	100/pk	R111672
	15x150mm - Double Bag	10/pk	R111675
Settling Plate Tryptic Soy Agar w/Lecithin, Polysorbate 80	15x150mm - Double Bag	10/pk	R111682

*Not intended for IVD use. [†] GripFit™ not available on these products.

©2009 Thermo Fisher Scientific. All rights reserved. All trademarks are property of Thermo Fisher Scientific and its subsidiaries.